SCALES

Digital Baby Scale

This digital baby scale is designed for accurate, dependable weighing with a 30 lb capacity. The scale features an easy-to-read LCD display, on-keypad instructions, and sleep and auto-shutoff modes to help extend battery life.



BB068

SPECIFICATIONS

Construction

- 6-digit, 7-segment 1" transflective LCD digital display
- On-keypad operating instructions
- Integrated infant weighing tray

Dimensions

- Overall: 26" W × 14.2" D × 7.5" H
- Tray: 26" W × 12.2" D

Utilities / Power

- (6) C-size alkaline, Ni-Cad, or NiMH batteries — sold separately
- 100–240 VAC power
- 12 VDC / 1.2 A medical UL-approved AC adapter

Performance

- Capacity: 30 lb / 15 kg
- Resolution: 0.1 oz / 5 g
- Sleep and auto-shutoff modes

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI