SWORDFISH PATHOLOGY SAW BLADES

Section Pathology Saw Blade

This section pathology saw blade is designed for both skullcap removal and smaller bone specimens. The wide portion supports long, linear cuts, while the smaller portion is designed for removal of small bone specimens. A coarse tooth pattern helps minimize blockage during use.



BD142

SPECIFICATIONS

Construction
- Section blade design
- Arborless mounting design

Materials
- Stainless steel

Performance
- Wide section designed for skullcap removal and long, linear cuts
- Smaller section designed for removal of small bone specimens
- Coarse tooth pattern helps minimize blockage
- Recommended blade replacement every 3 months

Dimensions
- 64 mm diameter

Compatibility
- Compatible with Mopec Swordfish Pathology Saw

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI