CLEANERS AND DISINFECTANTS

Eazy Soak Instrument Cleaner

Eazy Soak is a fast-acting, concentrated instrument cleaner designed to loosen and remove dried-on organic matter from laboratory instruments. Its low-foaming, mild-pH formula is suitable for hard-to-clean instruments such as scissors, forceps, serrated cutting blades, and scalpel handles.



BE046

SPECIFICATIONS

Performance

- Super concentrated
- Fast acting
- Low foaming
- Mild pH
- Helps loosen and remove dried-on organic matter, including blood, tissue, and soil
- Environmentally safe
- Makes up to 32 gallons of cleaning solutions
- Do not mix with bleach, peroxides, or other detergents

Capacity

- 1 quart

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI