CLEANERS AND DISINFECTANTS

SaniPath Disinfecting Spray Cleaner

SaniPath Disinfecting Spray Cleaner is a multi-purpose decontaminating cleaner designed for daily use on pathology and morgue work surfaces. It cleans, disinfects, and deodorizes while leaving behind a mild lemon fragrance.



BE047

SPECIFICATIONS

Performance

- Cleans, disinfects, and deodorizes
- Tuberculocidal and virucidal
- Suitable for daily use and spill cleanup
- Designed for equipment, tables, coolers, cutting stations, instrumentation, floors, computers, and telephones
- Mild lemon fragrance
- EPA Reg. No. 1839-83-74952
- EPA Est. No. 5747-MI-01

Capacity

- 32 fl oz

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI