
CLEANERS AND DISINFECTANTS

Formaway Formalin Neutralizing Absorbent Pads

Formaway Formalin Neutralizing Absorbent Pads are designed to reduce formaldehyde vapors by absorbing and neutralizing formaldehyde, glutaraldehyde, and other aldehydes. They provide vapor control during tissue handling, storage, and small formalin spills.



SPECIFICATIONS

Performance

- Absorbs and neutralizes aldehydes
- Reduces formaldehyde vapors
- Suitable for tissue handling, storage, and small spills
- Converts absorbed aldehydes to a non-hazardous polymer
- Blue color improves specimen visibility

Construction

- Soft cutting surface

Variations

Model No.	Description	Dimensions	Quantity
BE092	Mini Absorbent Pads	3" X 4"	100 per pack
BE093	Small Absorbent Pads	8" X 11"	25 per pack
BE094	Medium Absorbent Pads	11" X 15"	25 per pack
BE095	Large Absorbent Pads	15" X 20"	25 per pack
BE096	Absorbent Roll	13" X 30'	1 roll

MOPEC GROUP

Mopec

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI