CLEANERS AND DISINFECTANTS

EZ-Zyme Instrument Cleaner Packets

EZ-Zyme is an enzyme-based instrument cleaner designed to dissolve tissue debris from hard-to-clean dissection instruments. Each concentrated packet is diluted in warm water to create a cleaning solution that reaches cracks and crevices without hand scrubbing.



BE100

SPECIFICATIONS

Performance

- Dissolves tissue debris
- Reaches cracks and crevices
- No hand scrubbing required
- Designed for hard-to-clean instruments
- Enzyme-based cleaner
- Not a disinfectant
- 1 packet makes 1 gallon of cleaning solution

Quantity

- 32 packets

MOPEC GROUP
Mopec
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI