PERSONAL PROTECTIVE EQUIPMENT (PPE)

No-Contact Infrared Thermometer

This no-contact infrared thermometer is designed for body and forehead temperature screening without physical contact. Its three-color backlit display helps identify normal, elevated, and high temperature readings at a glance.



BE195

SPECIFICATIONS

Performance

- No-contact body and forehead temperature measurement
- Accuracy: ±0.3°C / ±0.6°F
- Stores up to 34 readings
- Color indicates normal, slightly high, and serious temperatures
- Switchable °C / °F units

Construction

- Three-color alarm backlit
- Handheld thermometer design

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI