GLOVES AND GLOVE DISPENSERS

Black Nitrile Exam Gloves, 6 Mil

These black nitrile exam gloves are designed for durable hand protection during laboratory and post-mortem procedures. Textured fingertips provide a secure grip, while the 6 mil nitrile construction offers resistance to oils, grease, and chemicals. Gloves are powder-free and available in sizes Small through XL.



SPECIFICATIONS

Performance

- Textured fingertips provide secure grip
- Resistant to oils, grease, and chemicals
- Designed for strength and durability
- Powder-free

Quantity

- Sizes S-L: 10 boxes (100 gloves per box) / case
  Size XL: 10 boxes (90 gloves per box) / case

Dimensions

- 6 mil thickness

Materials

- Nitrile
- Not made with natural rubber latex

Variations

Model No.	Description	Size	Quantity
BE501-S	Black Nitrile Exam Gloves	Small	100/box
BE501-M	Black Nitrile Exam Gloves	Medium	100/box
BE501-L	Black Nitrile Exam Gloves	Large	100/box
BE501-XL	Black Nitrile Exam Gloves	X-Large	90/box

MOPEC GROUP

Mopec | MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI