PREFILLED FORMALIN CONTAINER

Small Prefilled Formalin Containers

These small prefilled formalin containers are designed for collecting, storing, and transporting tissue specimens. Each polypropylene container is half-filled with 10% neutral buffered formalin and features a polyethylene screw cap with a secure closure. Available in multiple capacities and sold by the case.



SPECIFICATIONS

Construction
- Secure screw-cap closure

Quantity
- BG502: 128/case
- BG504: 50/case
- BG506: 50/case
- BG509: 50/case
- BG512: 50/case

Performance
- Designed for specimen collection, storage, and transport
- Secure closure helps maintain the container seal
- Supplied prefilled for convenient use

Capacity
- Available in multiple capacities

Materials
- Half-filled with 10% neutral buffered formalin
- Polypropylene
- Polyethylene

Variations

Model No.	Description	Capacity	Quantity
BG502	Prefilled Formalin Container	20 mL	128 per case
BG504	Prefilled Formalin Container	40 mL	50 per case
BG506	Prefilled Formalin Container	60 mL	50 per case
BG509	Prefilled Formalin Container	90 mL	50 per case
BG512	Prefilled Formalin Container	120 mL	50 per case

MOPEC GROUP
Mopec
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI