PREFILLED FORMALIN CONTAINER

Prefilled Formalin Container, 0.6 Gallon

This prefilled formalin container is half-filled with 10% neutral buffered formalin and features a tamper- and moisture-resistant screw-top seal. An integrated handle supports easy transport.



BG505

SPECIFICATIONS

Construction
- Screw-top closure
- Integrated carrying handle

Materials
- High-density polyethylene
- Half-filled with 10% neutral buffered formalin

Performance
- Designed for specimen collection, storage, and transport
- Tamper and moisture resistant seal
- Integrated handle supports easy transport
- Supplied prefilled for convenient use

Capacity
- 0.6 gallon

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI