PREFILLED FORMALIN CONTAINER

Prefilled Formalin Containers

These prefilled specimen containers are designed for collecting, storing, and transporting tissue specimens. Each polypropylene container is half-filled with 10% neutral buffered formalin and features a heat-sealed screw-cap closure to help prevent leakage during initial delivery.



SPECIFICATIONS

Construction
- Polypropylene container
- Polyethylene screw cap
- Heat-sealed lid

Materials
- Half-filled with 10% neutral buffered formalin

Performance
- Designed for specimen collection, storage, and transport
- Heat seal helps prevent leakage during initial delivery
- Available in multiple capacities

Quantity
- BG518: 50/case
- BG524: 12/case
- BG548: 24/case
- BG500: 12/case

Variations

Model No.	Description	Capacity	Quantity
BG518	Prefilled Formalin Container	180 mL	50 per case
BG524	Prefilled Formalin Container	240 mL	12 per case
BG548	Prefilled Formalin Container	480 mL	24 per case
BG500	Prefilled Formalin Container	1,000 mL	12 per case

MOPEC GROUP
Mopec
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI