ALCOHOL AND REAGENTS

Reagent Alcohol, 1 Gallon

This reagent alcohol is a blend of ethyl alcohol, isopropanol, and methanol intended for use as a general-purpose laboratory reagent. Available in multiple concentrations and supplied in 1-gallon containers, four per case.



SPECIFICATIONS

Performance

- General-purpose laboratory reagent
- Available in multiple concentrations

Materials

- Ethyl alcohol
- Isopropanol
- Methanol

Capacity

- 1 gallon per container

Quantity

- 4 containers per case

Variations

Model No.	Description	Concentration
BG732	Reagent Alcohol	50%
BG733	Reagent Alcohol	70%
BG734	Reagent Alcohol	80%
BG735	Reagent Alcohol	95%
BG736	Reagent Alcohol	100%

MOPEC GROUP

Mopec MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI