
HISTOLOGY EQUIPMENT

Mopec Heated Forceps System

The Mopec Heated Forceps System is designed for controlled specimen handling while minimizing bench-top space. The system includes a control unit, power supply, and two heated forceps, with either serrated or smooth jaw options.



SPECIFICATIONS

Construction

- Control unit and power supply
- Includes 1 mm heated forceps (Red)
- Includes 2 mm heated forceps (Yellow)

Materials

- Stainless steel forceps
- Silicone finger grips

Performance

- Operates at 145–158°F (63–70°C)
- Supports controlled specimen placement
- Compact bench-top design
- Control unit supports up to two heated forceps or heated tampers
- Available with serrated or smooth jaws

Dimensions

- Control unit: 120 × 120 × 40 mm
- Forceps length: 150 mm
- Tip length: 19.5 mm
- Cable length: 1.1 m

Utilities / Power

- 100/240 V AC
- Interchangeable US, UK, and EU sockets

Variations

Model No.	Description	Jaw Style
BK300-SE	Heated Forceps System	Serrated
BK300-SM	Heated Forceps System	Smooth

MOPEC GROUP

Mopec

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI