MICROTOMY AND CYTOLOGY

Mopec Microtome Blades

These microtome blades feature a composite coating designed to produce high-quality tissue sections while supporting extended blade life. Each blade is tested for compression, cutting feel, and longevity and is available in low- and high-profile configurations.



BM100-LP



BM100-HP

SPECIFICATIONS

Performance

- Designed for high-quality tissue sectioning
- Composite coating supports extended blade life
- Tested for compression, cutting feel, and longevity

Construction

- Disposable microtome blade
- Available in low- and high-profile configurations
- Composite-coated blade

Quantity

- 50 blades per pack

Dimensions

- BM100-LP: Low Profile: 3″ L x .315″ W x .010″ Thick
- BM100-HP: High Profile: 3″ L x .554″ W x .017″ Thick

Variations

Model No.	Description	Profile
BM100-LP	Microtome Blades	Low Profile
BM100-HP	Microtome Blades	High Profile

MOPEC GROUP

Mopec

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI