GROSSING STATION ACCESSORIES

Maestro Hand Spray Rinse



MB1032R



MB1032L

SPECIFICATIONS

Construction

- Designed for customers to rinse off the workstation of their grossing station
- Includes (1) sprayer, (1) 60" hose, flow control fixture, plumbing and fittings
- Placed on the opposite side of the unit's sink
- Now standard in all Maestro Encore grossing stations

Dimensions

- 60" hose

Materials

- Hand sprayer with soft spray: chrome-plated brass
- Sheetmetal: 304 stainless steel construction with a polished No. 4 finish
- Hardware: 304 stainless steel
- Plumbing lines and fittings: brass, pex, and stainless steel

Performance

- Mixing valve/flow control fixture
- Backflow preventers on both hot and cold-water lines

Options & Accessories

- MB1032L: designed for left-hand Maestro grossing stations
- MB1032R: designed for right-hand Maestro grossing stations

Compatibility

- Compatible with all Maestro grossing stations that do not currently have this option

NOTE:

- For stations adding this option where it had not previously existed, it would include a new endcap to the worksurface with a hole for the hand spray to retract inside the station.

Variations

Model No.	Description	Compatibility
MB1032L	Hand Spray Rinse, Left-Hand Units	Left-Hand Maestro Grossing Stations
MB1032R	Hand Spray Rinse, Right-Hand Units	Right-Hand Maestro Grossing Stations

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI