GROSSING STATION ACCESSORIES

Ensemble Pegboard Barcode Holster



SPECIFICATIONS

Construction
- Designed for the Mopec Maestro's Ensemble Pegboard
- Designed for the most popular barcode reader guns

Capacity
- 4 lbs

Materials
- 304 stainless steel construction with a polished No. 4 finish and ABS plastic

Dimensions
- 6" L X 7.4" D X 13.6" H

Compatibility
- Fits all Mopec Pegboards, including AX031R and AX031L pegboard wings
- Can also be attached to the front of grossing stations, above the vent grille

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI