Instrument Tray & Lid Option With Flush Handle

This instrument tray features a flush handle and optional lid for convenient storage and soaking of soiled dissection instruments during laboratory procedures. Its practical design allows for easy handling and organization, providing a reliable solution for everyday laboratory use.



BA033



BA035



BA034



BA039

SPECIFICATIONS

Materials

Stainless steel

Dimensions

- 20.75" (528 mm) L X 12.75" (324 mm) W
- BA033: 2.5" D
- BA034: 4" D
- BA035: 6" D

Compatibility

- BA039: Universal Cover is made to fit all trays
- For use with Instrument Cleaner (BE046)
- NOTE: Unless choosing the Tray & Lid Bundle, the items are all sold separately.

Variations

Model No.	Description	Depth
BA033	Instrument Tray	2.5" D
BA034	Instrument Tray	4" D
BA035	Instrument Tray	6" D
BA039	Instrument Tray Cover	Lid

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com

800.362.8491

800 Tech Row, Madison Heights, MI