Large Specimen Bowl Strainer

This strainer allows for slow, controlled straining and rinsing. It is made of stainless steel for durability and easy cleaning. Its large size provides space for handling specimens during processing. It is designed for use in pathology and laboratory settings.



BA050

SPECIFICATIONS

Construction

- Approximately 300 holes, 0.25″ in diameter

Materials

- Stainless steel

Dimensions

- 16″ (406 mm) Diameter X 7″ (178 mm) Deep

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI