BUCKETS, TRAYS AND STRAINERS

Instrument & Specimen Utility Bin With Strainer

This perforated multipurpose bin with strainer is designed for temporary specimen fixation, instrument soaking, and specimen collection during laboratory procedures. Its side handles allow for easy lifting, efficient drainage, and rinsing, while its stackable design provides convenient storage.



BA053

SPECIFICATIONS

Construction

- Approximately 300 holes, 0.25″ diameter

Materials

- Polyethylene

Dimensions

- Top: 20.5″ (52.07 cm) L x 16″ (40.64 cm) W
- Strainer: 20″ (50.80 cm) L x 15″ (38.1 cm) W x 5″ (12.7 cm) D
- Base: 20″ (50.80 cm) L x 15.5″ (39.37 cm) W x 7″ (17.78 cm) D

Compatibility

- Can be used with Instrument Cleaner (BE046)

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI