



REV	DESCRIPTION	DATE	ECN #	APPR
A	INITIAL RELEASE	00-XXX-19	XXXXX	MM





DETAIL A
(4 PLACES)





Item #	Qty.	File Name	Title	Material	Rev.
1	4	pb0013.par	NUT - #6-32, NYLON LOCK, 18-8 SS	PURCHASED PART (SEE PART DESCRIPTION)	A
2	4	pb0027.par	SCREW - SS, 6-32 x 1/2"L, 82 DEG FLAT PHILLIPS	PURCHASED PART (SEE PART DESCRIPTION)	A
3	4	pm0166.par	PM0166.par	PURCHASED PART (SEE PART DESCRIPTION)	A
4	1	S2159-A.psm	NECK BRACE S.S.	SHT, SS, 16 GA, 304	A

MOPEC LEGAL NOTE
THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION OF MOPEC, AND MAY NOT BE COPIED IN WHOLE OR PART, AND NO USE SHALL BE MADE OF IT OTHER THAN FOR WHICH IT HAS BEEN SUBMITTED WITHOUT THE WRITTEN PERMISSION OF A MOPEC MANAGER
MOPEC 2026

UNLESS OTHERWISE SPECIFIED:
ALL DIMENSIONS IN INCHES
TOLERANCES
0.0 ± 0.1 in
0.00 ± 0.05 in
0.000 ± 0.010 in
0.0000 ± 0.0010 in
ANGLES ±0.5°
BREAK ALL CORNERS 0.063 X 45°
ALL FILLET RADII R0.063 /0.125in
INSPECT OVAL DIMENSIONS

MATERIAL:

MOPEC GROUP
Mopec MORTECH MANUFACTURING
DRAWN TO 3RD ANGLE PROJECTION

DRAWN BY: D. ANDERSON
DATE: 11/3/26
CHECKED BY: bdevins
DATE: 11/3/2026
APPROVED BY:
DATE: ---

PART NUMBER: BA078-A
DESCRIPTION: NECK REST
SIZE: A2 FILE NAME: BA078-A.dft REV: A
SCALE: 1:2
DO NOT SCALE IF IN DOUBT ASK
SHEET 1 OF 1

BA078-A.dft