SCALES

Digital Compact Scale

This compact digital scale features a lightweight, durable, and portable design for convenient weighing in a variety of laboratory applications. Its broad weighing platform accommodates larger specimens, providing reliable measurements and versatility during routine laboratory procedures.



SPECIFICATIONS

Construction
- User calibration from keypad
- 11 mm LCD display

Dimensions
- Overall: 5.3" W X 7.9" D X 1.6" H
- Pan: 5" W X 5.6" D

Capacity
- BB061: 2,000 g
- BB062: 5,000 g

Compatibility
- (BB203) AC adapter

Utilities / Power
- Requires (3) AA alkaline batteries (included) or AC adapter (sold separately)

Performance
- Weighing units: g, lb, oz, newton

Variations

Model No.	Description	Capacity	Dimensions
HB1100	Digital Compact Scale	2,000 g	5.3" W X 7.9" D X 1.6" H
HB1410	Digital Compact Scale	5,000 g	5.3" W X 7.9" D X 1.6" H

MOPEC GROUP
Mopec ELEVATING PATHOLOGY
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI