CLEANERS AND DISINFECTANTS

ClearSteel Stainless Steel Cleaning & Polishing Wipes

Mopec's ClearSteel Wipes effortlessly clean and polish all stainless steel laboratory equipment. With the swipe of a towelette, stains, smudges and blood are eliminated. Surfaces will be left with a natural luster and a protective covering that repels water, fingerprints and dust. Specifically designed to be used on all equipment such as workstations, cabinets, coolers and storage areas.



BE039

SPECIFICATIONS

Construction
- Pop-up canister

Dimensions
- 6" L X 6" W

Performance
- Non-abrasive
- Repels water and dust

Quantity
- (160) premoistened towelettes

Note:
- **This product is not available in California.**

MOPEC GROUP
Mopec ELEVATING PATHOLOGY
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI