
GLOVES AND GLOVE DISPENSERS

Blue Nitrile Exam Gloves, 5 Mil

These blue nitrile exam gloves are designed for medical-grade hand protection during laboratory and postmortem procedures. Their 5 mil nitrile construction offers enhanced puncture resistance and flexibility, while the fully textured surface provides a secure grip in both wet and dry conditions.



SPECIFICATIONS

Construction
- Blue color

Materials
- Nitrile

Dimensions
- 5 mil thickness

Quantity
- Sizes XS-L: 10 boxes (100 gloves per box) / case
- Sizes XL-XXL: 10 boxes (90 gloves per box) / case

Performance
- Textured fingertips provide secure grip in dry or wet conditions
- Powder-free
- Improved puncture resistance

Variations

Model No.	Description	Size	Quantity
BE500-XS	Blue Nitrile Exam Gloves	X-Small	100/Box
BE500-S	Blue Nitrile Exam Gloves	Small	100/Box
BE500-M	Blue Nitrile Exam Gloves	Medium	100/Box
BE500-L	Blue Nitrile Exam Gloves	Large	100/Box
BE500-XL	Blue Nitrile Exam Gloves	X-Large	90/Box
BE500-XXL	Blue Nitrile Exam Gloves	XX-Large	90/Box

MOPEC GROUP
Mopec ELEVATING PATHOLOGY
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI