TISSUE MARKING SPRAYS

Amplify Setting Solution

This Amplify Setting Solution is a fixative spray designed to improve the adhesion of Harmony Tissue Marking Dyes on challenging specimens. Its specialized formula modifies tissue surfaces and can be applied before or after dye application.



BG100

SPECIFICATIONS

Performance

- Modifies the surface of challenging specimens
- Aids the adhesion of Harmony Tissue Marking Dyes
- Can be applied to tissue before its marked with dye
- Suitable for fatty specimens

Capacity

- 8 oz (250 ml) bottle

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI