FORMALIN NEUTRALIZERS

Formaway Green Liquid Formalin Neutralizer

Mopec's Formaway Green Liquid Formalin Neutralizer quickly and effectively neutralizes spent formalin for sewer disposal. Its convenient liquid formula provides a practical solution for managing formalin waste in pathology and laboratory environments. Available in 1 gallon, 2.5 gallon, and 5 gallon cubes.



SPECIFICATIONS

Performance

- Add 32 ounces of Formaway Green to each gallon of 10% neutral buffered formalin waste and allow to stand for at least 4 hours before disposing

Capacity

- BG560-1: 1 gallon cube
- BG560-2.5: 2.5 gallon cube
- BG560-5: 5 gallon cube

Variations

Model No.	Description	Capacity
BG560-1	Formaway Green Formalin Neutralizer	1 gallon
BG560-2.5	Formaway Green Formalin Neutralizer	2.5 gallon
BG560-5	Formaway Green Formalin Neutralizer	5 gallon

MOPEC GROUP
Mopec ELEVATING PATHOLOGY
MORTECH MANUFACTURING
www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI