TISSUE FLOTATION WATER BATH

Mopec Illuminated Tissue Flotation Water Bath

Mopec's Illuminated Tissue Flotation Water Bath features an LED light and high-contrast background that allows users to better view tissue sections while also detecting foreign particles floating on the surface of the water. The LED light is dimmer controlled and can be switched between white or blue.



BK200



BK200-L

SPECIFICATIONS

Construction

- Switch between blue and white LED back lighting with dimmer control for optimal viewing
- Extra deep bowl to assist with tissue pick-up at any angle
- Digital temperature control
- Convenient slide drying area at the rear of the bath
- Anti-microbial coating to help prevent bacteria and virus growth

Dimensions

- Overall: 295 mm W X 370 mm D X 135 mm H
- Bowl: 90 mm H X 220 mm Diameter

Utilities / Power

- Dual Voltage; Operates under both 110V and 230V power supplies

Compatibility

- Optional lid helps improve temperature control and prevent surface particles (BK200-L)

Variations

Model No.	Description
BK200	Illuminated Tissue Flotation Bath
BK200-L	Illuminated Tissue Flotation Bath Lid

MOPEC GROUP

Mopec ELEVATING PATHOLOGY

MORTECH MANUFACTURING

www.mopecgroup.com
800.362.8491
800 Tech Row, Madison Heights, MI